




NEWCLIP
TECHNICS



SURGICAL TECHNIQUE

Hyla
HYBRID LATARJET


SURGICAL TECHNIQUE

Example of a surgical procedure



1. Perform the release and the osteotomy of the coracoid process.

Prepare the coracoid process, then measure its width using the double ruler (ANC919).



2. Directly report the width measured to set the intra-articular bracket offset (ANC913) using the screwdriver (ANC917) with the quick coupling handle (ANC350).

This step allows an optimized centering of the screws into the coracoid process.



3. Temporarily stabilize the coracoid process onto the guide using the holding instrument (ANC914) by maintaining traction with the hand (see highlighted area).

For more information on how to assemble the instrument, see page 3.




4. While holding the coracoid process, slide a 2-in-1 instrument (ANC915) in each cannulas of the guide until they reach the coracoid process.



5. Slide a drill bit (ANC918) into the cannula of one of the 2-in-1 instruments (ANC915) (a) and insert a Ø1.2 mm pin (33.0212.200) (b) through its cannula and into the coracoid process.

Then drill the first cortex using the drill bit (ANC918) still in position (c).



6. Using the cannulated handle (ANC143), with the drill bit and pin still in place, completely screw the 2-in-1 instrument (ANC915) into the coracoid process to achieve prehension (a).

Remove the drill bit and the pin.

Repeat steps 5 and 6 for the other 2-in-1 instrument (ANC915) (b).



7. Once both 2-in-1 instruments are completely inserted make sure they are in contact with the rear of the guide - highlighted area). Release and remove the holding instrument by sliding it backwards.


Check that the coracoid process is attached to the instrument.

SURGICAL TECHNIQUE





8. Position the guide-coracoid process assembly against the anterior part of the glenoid rim. The intra-articular bracket (ANC913) must be positioned in contact with the articular glenoid surface.

The intra-articular bracket allows to ensure the appropriate medio-lateral positioning of the coracoid process.



9. Once the appropriate position is obtained, slide a drill bit (ANC918) into the cannula of the inferior 2-in-1 instrument and insert a 1.2 mm pin (33.0212.200) until it reaches the posterior cortex of the scapula. Following the same procedure, insert another pin in the superior cannula.

At the rear of the guide, position the ruler in order to determine the screw length required. Directly read the markings at the rear of the pins. The ruler can be positioned in a vertical or horizontal position to ease the reading.





10. Through the inferior 2-in-1 instrument (ANC915), drill the first cortex of the scapula using the drill bit (ANC918) still in place (a).

Remove the drill bit. Remove the 2-in-1 instrument (ANC915) using the cannulated handle (ANC143) and insert a Ø3.5 mm cannulated screw using the screwdriver (ANC917) with the quick coupling handle (ANC350). Remove the pin (b).




11. Repeat the step 10 for the superior 2-in-1 instrument.

FINAL RESULT



REFERENCES

INSTRUMENTATION		
Ref.	Description	Qty
ANC350	Ø4.5 mm AO quick coupling handle - Size 1	2
ANC912	Double targeting device for Latarjet technique	1
ANC913	Offset intra-articular bracket	1
ANC914	Holding instrument	1
ANC915	3-in-1 instrument for Latarjet technique - Short	1
ANC916	3-in-1 instrument for Latarjet technique - Long	1
ANC917	2.5 mm quick coupling hexagonal non prehensor screwdriver - cannula Ø1.3 mm	1
ANC918	Ø2.5 mm quick coupling drill bit - cannula Ø1.3 mm	1
ANC919	Double ruler	1
33.0212.200	Pin Ø1.2 L200 mm	4



IMPLANTS	
Ref.	Description
H1.25IFT3.5Lxx-ST	Self-drilling self-compressive screw Ø3.5 mm - cannulated Ø1.25 mm - STERILE L28 mm to L42 mm (2 mm increment)



S-Box : sterile screws

The information presented in this brochure is intended to demonstrate a NEWCLIP TECHNICS product. Always refer to the package insert, product label and/or user instructions before using any NEWCLIP TECHNICS product. Surgeons must always rely on their own clinical judgment when deciding which products and techniques to use with their patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your NEWCLIP TECHNICS representative if you have questions about the availability of NEWCLIP TECHNICS products in your area.

NEWCLIP TECHNICS

PA de la Lande Saint Martin
45 rue des Garottières
44115 Haute Goulaine, France
+33 (0)2 28 21 23 25
orders@newcliptechnics.com
www.newcliptechnics.com

NEWCLIP TECHNICS GERMANY

Newclip GmbH
Pröllstraße 11
D-86157 Augsburg, Germany
+49 (0)821 650 749 40
info@newclipgmbh.com
www.newclipgmbh.de

NEWCLIP TECHNICS USA

Newclip USA
340 Tesconi Circle Suite A
Santa Rosa CA 95401, USA
+1 707 230 5078
customerservice@newclipusa.com
www.newclipusa.com

NEWCLIP TECHNICS AUSTRALIA

Newclip Australia
3B/11 Donkin Street
West End 4101, Australia
+61 (0)2 81 886 110
solutions@newclipaustralia.com
www.newcliptechnics.com

NEWCLIP TECHNICS JAPAN

Newclip Technics Japan K.K.
KKK Bldg. 502, 3-18-1 Asakusabashi
Taito-Ku, Tokyo, 111-0053, Japan
+81 (0)3 58 25 49 81
Fax: +81 (0)3 58 25 49 86
www.newcliptechnics.com

NEWCLIP TECHNICS IBERIA

Newclip Iberia
Calle Frederic Mompou, 4b
Sant Just Desvern, 08960 Barcelona, Spain
+34 938 299 526
contact@newclipiberia.com
www.newcliptechnics.com



NEWCLIP-TECHNICS

MADE IN FRANCE