

SINGLE USE KIT

STERILE R



TM

NEWCLIP  
TECHNICS



INITIAL

Ankle



Ready  
when you are!

With a non sterile standard kit



Calling on medical staff

Constraints



Complex traceability



Contracted out sterilization




Suppliers' deadline


### High costs



- \$ Stocks
- \$ Control
- \$ Cleaning
- \$ Decontamination
- \$ Sterilization



Bulky storage



### Complex process



Prevents an effective solution & a quick response



Defective sterilization



Incomplete kit



Damaged instrumentation



**INCREASED RISKS**  
**NON OPTIMIZED surgery**



**URGENT SURGICAL CASES COMPROMISED**



Safety



TRACEABILITY  
**100%**



Risk of contamination


### Cost efficiency

Controlled stocks  
Simplified control

- 0 Cleaning
- 0 Decontamination
- 0 Sterilization

Sundry expenses

Optimized storage



### Efficiency

1 Delivery

2 Storage

3 Surgery

An effective solution & a quick response

### Available when needed

READY-TO-USE FOR SURGERY

Optimized handling of URGENT SURGICAL CASES

STERILE R SINGLE USE KIT  
with state-of-the-art implants

Ready when you are!



### Safety:

The Initial A™ kits are fully traceable and have a shelf life of 5 years. Its instrumentation and implants are “always new” and have never been opened or used before.



### Available when needed:

The Initial A™ kits (Initial A™ - Fibula and Initial A™- Syndesmosis) come pre-sterilized and ready to use. The combination of sterile implants and single use instrumentation in a single packaging makes Initial A™ ideal for use in urgent surgical cases.



### Storage:

Initial A™ kits can be easily stored in the operating room because of its small size.



### Costs:

Initial A™ is a cost-effective solution. The additional costs including cleaning, decontamination, sterilization of kits are cancelled.



### Contamination:

The combination of sterile implants and sterile single-use instrumentation minimizes contamination risks.



### Buying procedure:

Initial A™ facilitates buying procedures: restocking and orders are simplified, stock management is optimized.

# Kit Content




## > Indications

The implants of the Initial A™ range are intended for the fixation of fractures, osteotomies and pseudarthroses of the distal and the diaphyseal fibula, the distal tibia and for the syndesmotomic repair in adults.


## > Contraindications

- Serious vascular deterioration, bone devitalization.
- Pregnancy.
- Acute or chronic, local or systemic infections.
- Lack of musculo-cutaneous cover, severe vascular deficiency affecting the concerned area.
- Insufficient bone quality preventing a good fixation of the implant into the bone.
- Muscular deficit, neurological deficiency or behavioural disorders which could submit the osteosynthesis to abnormal mechanical strains.
- Allergy to one of the materials used or sensitivity to foreign bodies.
- Serious problems of non-compliance, mental or neurological disorders, failure to follow post-operative care recommendations.
- Unstable physical and/or mental condition.

## > Initial A™ - Fibula kit



## > Initial A™ - Syndesmosis kit




# Plate features

## > A comprehensive range of plates




### Standard plates

Fixation of osteoporotic bones and complex fractures with or without syndesmosis injuries (green anodized plates for the right side, blue anodized plates for the left side).





### Narrow plates


Fixation of simple fractures with or without syndesmosis injuries (same plate for both sides).



## > Fixation system



Locking screw - Ø3.5 mm  
Ref. SOT3.5Lxx



Standard cortical screw - Ø3.5 mm  
Ref. CT3.5Lxx



Locking screw - Ø2.8 mm  
Ref. SDT2.8Lxx




Lag screw - Ø4.0 mm  
Ref. QT4.0Lxx



# Technical features: Initial A - Fibula kits

## > Precontoured implant


The design of this implant is the result of a proprietary state-of-the-art mapping technology to establish an optimized congruence between the plate and the bone.



### ► PLATE BENDING

The implant also offers bendable areas which allows adjusting of the plate on the diaphyseal part and on the junction of the diaphysis and epiphysis parts thanks to the bending pliers. **They are available separately, on demand, in non sterile version.**

Bending is only possible in the areas intended for this purpose. A bendable area must be bent only once, in one direction and not be performed excessively. The holes must be protected so as to avoid damaging the fixation.




**Handle for guide gauge:** before performing the drilling into the oblong hole, clip the handle for guide gauge on the Ø2.7 mm threaded guide gauge.



## > Angular range: +/- 10° polyaxial locking fixation


Initial A™ plates combine both polyaxial and locking technologies to create a fixed-angle construct particularly useful for poor bone quality and/or multifragmentary fractures.



Dualtec System® II  
Technology  
Clip + nut

## > Monoaxial locking system

- The threaded sections under the screw head and inside the hole have the same characteristics (1),
- Screw head cap (2),
- Implants material: titanium alloy.




# Surgical technique: Initial A - Fibula kits

Example: surgical technique with a standard size 2 right (KIT-AL2D).

(Same technique for all standard and narrow plates)


Page 1/2

## 2-IN-1 INSTRUMENT (ANC543)




The 2-in-1 instrument includes the 2 following functions:

- **Screwdriver** for  $\varnothing 3.5$  mm and  $\varnothing 4.0$  mm screws,
- **Countersink** to widen the drilling made in the first cortex before screw insertion.



**1.** Using the template (ANC607 or ANC659), define the suitable plate size, then determine the appropriate kit.


*N.B.: The templates can be used both for the right side and for the left side and are available separately in a sterile version.*




**2.** To insert an interfragmentary screw, drill using the  $\varnothing 2.7$  mm drill bit.



**3.** When a lag effect is necessary, use the countersink part of the blue 2-in-1 instrument to widen the first cortex previously drilled.




**4.** Insert the interfragmentary light blue  $\varnothing 3.5$  mm cortical screw using the screwdriver part of the blue 2-in-1 instrument.



**5.** Hold the plate by inserting pins through the dedicated distal holes.

The pins can be removed once the plate is stabilized.



**6.** Clip the handle for guide gauge on the blue  $\varnothing 2.7$  mm threaded guide gauge and perform the drilling using the assembly in the most distal oblong hole.

**Option 1** - Determine the screw length using the drill bit and guide gauge.


**Option 2** - Determine the screw length using the length gauge.

# Surgical technique: Initial A - Fibula kits


Example: surgical technique with a standard size 2 right (KIT-AL2D).

(Same technique for all standard and narrow plates)

Page 2/2




**7.** Screw a light blue Ø3.5 mm cortical screw in the oblong hole using the screwdriver part of the blue 2-in-1 instrument to secure the plate in place.




**8.** Using the grey Ø2.0 mm threaded guide gauge, choose the angle of the non-anodized Ø2.8 mm locking screws in the polyaxial holes then drill (Ø2.0 mm).

**Option 1** - Determine the screw length using the drill bit and guide gauge.

**Option 2** - Determine the screw length using the length gauge.




**9.** Using the grey screwdriver, insert and lock the non-anodized Ø2.8 mm locking screws.




**10.** Using the blue Ø2.7 mm threaded guide gauge, drill (Ø2.7 mm).

**Option 1** - Determine the screw length using the drill bit and guide gauge.

**Option 2** - Determine the screw length using the length gauge.



**11.** Using the countersink part of the blue 2-in-1 instrument, widen the first cortex previously drilled. Insert a blue Ø3.5 mm locking screw using the screwdriver part of the blue 2-in-1 instrument and lock it.



Repeat previous steps to insert the remaining Ø3.5 mm screws in the plate.

**⚠** The final tightening of the screws must be performed by hand.

# References: Initial A - Fibula kits



INITIAL A™ - FIBULA KITS	
Ref.	Description
KIT-AL1D	Distal fibula kit - Standard - Right - Size 1
KIT-AL1G	Distal fibula kit - Standard - Left - Size 1
KIT-AL2D	Distal fibula kit - Standard - Right - Size 2
KIT-AL2G	Distal fibula kit - Standard - Left - Size 2
KIT-AL3D	Distal fibula kit - Standard - Right - Size 3
KIT-AL3G	Distal fibula kit - Standard - Left - Size 3
KIT-AL1S	Distal fibula kit - Narrow symmetrical - Size 1
KIT-AL2S	Distal fibula kit - Narrow symmetrical - Size 2

INITIAL A™ - FIBULA KIT - INSTRUMENTATION CONTENT	
Description	
Ø2.0 mm quick coupling drill bit - L 125 mm	
Ø2.7 mm quick coupling drill bit - L 125 mm	
Ø2.0 mm threaded guide gauge for Ø2.8 mm screws	
Ø2.7 mm threaded guide gauge for Ø3.5 mm screws	
Handle for guide gauge	
Length gauge for Ø2.8 and Ø3.5 mm screws - L 10-60 mm	
2.0 mm hexagonal prehensor screwdriver	
2 in 1 : 2.5 mm hexagonal prehensor screwdriver - Ø3.5 mm countersink	
Pin - Ø1.4 L120 mm (x2)	


NB: Supplemental screws are available in sterile package (cf: Initial A™ additional kits, additional implants).


INITIAL A™ - FIBULA KITS - IMPLANTS CONTENT			QUANTITY PER KIT				
	Ref.	Description	KIT-AL1D or KIT-AL1G	KIT-AL2D or KIT-AL2G	KIT-AL3D or KIT-AL3G	KIT-AL1S	KIT-AL2S
STANDARD PLATES	RTDLS1 or RTGLS1	Lateral plate for distal fibula - Standard Right or left - Size 1	1	-	-	-	-
	RTDLS2 or RTGLS2	Lateral plate for distal fibula - Standard Right or left - Size 2	-	1	-	-	-
	RTDLS3 or RTGLS3	Lateral plate for distal fibula - Standard Right or left - Size 3	-	-	1	-	-
NARROW PLATES	RTSLN1	Lateral plate for distal fibula - Narrow symmetrical - Size 1	-	-	-	1	-
	RTSLN2	Lateral plate for distal fibula - Narrow symmetrical - Size 2	-	-	-	-	1
LOCKING SCREWS Ø2.8 MM	SDT2.8L10	Locking screw - Ø2.8 mm - L 10 mm	1	1	1	-	-
	SDT2.8L12	Locking screw - Ø2.8 mm - L 12 mm	1	1	1	-	-
	SDT2.8L14	Locking screw - Ø2.8 mm - L 14 mm	2	2	2	1	1
	SDT2.8L16	Locking screw - Ø2.8 mm - L 16 mm	2	2	2	2	2
	SDT2.8L18	Locking screw - Ø2.8 mm - L 18 mm	2	2	2	1	1
LOCKING SCREWS Ø3.5 MM	SOT3.5L12	Locking screw - Ø3.5 mm - L 12 mm	1	2	3	2	2
	SOT3.5L14	Locking screw - Ø3.5 mm - L 14 mm	1	2	2	2	2
	SOT3.5L16	Locking screw - Ø3.5 mm - L 16 mm	1	1	1	1	2
	SOT3.5L18	Locking screw - Ø3.5 mm - L 18 mm	1	1	1	-	-
STANDARD CORTICAL SCREWS Ø3.5 MM	CT3.5L12	Standard cortical screw - Ø3.5 mm - L 12 mm	-	1	1	-	1
	CT3.5L14	Standard cortical screw - Ø3.5 mm - L 14 mm	1	1	2	1	1
	CT3.5L16	Standard cortical screw - Ø3.5 mm - L 16 mm	1	1	1	1	1
	CT3.5L18	Standard cortical screw - Ø3.5 mm - L 18 mm	-	1	1	-	-
	CT3.5L20	Standard cortical screw - Ø3.5 mm - L 20 mm	1	-	-	-	-
	CT3.5L22	Standard cortical screw - Ø3.5 mm - L 22 mm	-	1	1	1	1
	CT3.5L24	Standard cortical screw - Ø3.5 mm - L 24 mm	1	-	-	-	-

# Surgical technique: Initial A - Syndesmosis kit

Example: surgical technique with a right standard plate, size 2 (KIT-AL2D + KIT-AMS).


(Same technique for all standard and narrow plates)

 The final tightening of the screws must be performed by hand.




1. The syndesmosic screw can be inserted in the following holes:

- For the standard plates:
  - a. The most distal diaphyseal hole,
  - b. The most distal oblong hole,
- For the narrow plates:
  - c. The most distal oblong hole.



2. Drill using the non threaded bent long guide gauge for syndesmosic screws. The drilling length can be directly measured on the guide gauge.

 It is compulsory to use this guide.



**Final Result**

**Insert the syndesmosic screw and finalize the tightening with the screwdriver part of the blue 2-in-1 instrument from the Initial A™ - Fibula kit.**



*NB: The syndesmosic screw must be removed using the removal kit for Ø3.5 mm screws (ref: KIT-REMOVE-A) once the syndesmosis has healed, usually after six to eight weeks.*

# References: Initial A - Syndesmosis kit



## INITIAL A™ - SYNDESMOSIS KIT

Ref.	Description
KIT-AMS	Syndesmosis kit

## INITIAL A™ - SYNDESMOSIS KIT CONTENT


Description	Qty
Ø2.7 mm quick coupling drill bit - L 180 mm	1
Ø2.7 mm non threaded bent long guide gauge for Ø3.5 and Ø4.0 mm screws	1
2.5 mm quick coupling hexagonal screwdriver	1
Washer	2

## STERILE SCREWS FOR INITIAL A™ - AVAILABLE SEPARATELY - SYNDESMOSIS KIT\*


Ref.	Description	Qty
CT3.5L40-ST	Standard cortical screw - Ø3.5 mm - L 40 mm - STERILE	2
CT3.5L45-ST	Standard cortical screw - Ø3.5 mm - L 45 mm - STERILE	3
CT3.5L50-ST	Standard cortical screw - Ø3.5 mm - L 50 mm - STERILE	3
CT3.5L55-ST	Standard cortical screw - Ø3.5 mm - L 55 mm - STERILE	3
CT3.5L60-ST	Standard cortical screw - Ø3.5 mm - L 60 mm - STERILE	2
CT3.5L65-ST	Standard cortical screw - Ø3.5 mm - L 65 mm - STERILE	2
QT4.0L40-ST	Lag screw - Ø4.0 mm - L 40 mm - STERILE	3
QT4.0L45-ST	Lag screw - Ø4.0 mm - L 45 mm - STERILE	3
QT4.0L50-ST	Lag screw - Ø4.0 mm - L 50 mm - STERILE	3
QT4.0L55-ST	Lag screw - Ø4.0 mm - L 55 mm - STERILE	3
QT4.0L60-ST	Lag screw - Ø4.0 mm - L 60 mm - STERILE	3
QT4.0L65-ST	Lag screw - Ø4.0 mm - L 65 mm - STERILE	3

\*CT3.5Lxx: Light blue anodized.  
QT4.0Lxx: Non anodized.

# Surgical technique Initial S™ - 4.0 cannulated screw for medial malleolus fractures



1. Insert the Ø1.3 mm pin to stabilize the two fragments.




2. Slide the length gauge along the Ø1.3 mm pin until the cortex is reached (a). Determine the insertion depth using the marking on the pin (b).

NB: The pin can then be inserted deeper in order to prevent its removal during drilling.



3. Select the appropriate screw length and insert the screw along the pin using the screwdriver part of the 2-in-1 instrument until the desired reduction and compression are achieved. Then remove the pin.



**FINAL RESULT**



## Compression washer




NB: In case of osteoporotic bone, it is possible to add a compression washer under the screw head before step 3 to obtain an optimized compression.



## Optional steps:

These steps can be done before screwing.

1- In case of a hard bone density or several cortices, it is recommended to drill before the screw insertion. The drilling depth can be checked using the marking on the drill bit.



2- If reaming is required, widen the surface of the insertion using the countersink part of the 2-in-1 instrument.




## References

STERILE INSTRUMENTATION FOR Ø4.0 MM CANNULATED SCREWS			
Kit	Description	Content	Qty
KIT-SCQ4.0	Kit for Ø4.0 mm cannulated headed screws	Length gauge for pin Ø1.3 mm - L120 mm	1
		2 in 1: 2.5 mm hexagonal screwdriver - Ø6.0 mm countersink	1
		Ø2.9 mm drill bit - cannula 1.4 mm - L 120 mm - AO Ø4.5 mm quick coupling	1
		Washer	2
		Pin Ø1.3 L140 mm	3
		5.8 mm single use handle	1

CANNULATED HEADED SCREWS Ø4.0 MM*	
Réf.	Description
H1.4QT4.0Lxx-ST	Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - short thread - L xx mm (from 26 to 70) - 2 mm increment - STERILE

\*Available separately in sterile version.



# References: Additional kits


## > Additional implants

Sterile screws




LOCKING SCREWS - Ø2.8 mm*			
Ref.	Description	Qty	
SDT2.8L10-ST	Locking screw - Ø2.8 mm - L 10 mm - STERILE	1	
SDT2.8L12-ST	Locking screw - Ø2.8 mm - L 12 mm - STERILE	1	
SDT2.8L14-ST	Locking screw - Ø2.8 mm - L 14 mm - STERILE	1	
SDT2.8L16-ST	Locking screw - Ø2.8 mm - L 16 mm - STERILE	1	
SDT2.8L18-ST	Locking screw - Ø2.8 mm - L 18 mm - STERILE	1	
SDT2.8L20-ST	Locking screw - Ø2.8 mm - L 20 mm - STERILE	2	
SDT2.8L22-ST	Locking screw - Ø2.8 mm - L 22 mm - STERILE	2	
SDT2.8L24-ST	Locking screw - Ø2.8 mm - L 24 mm - STERILE	1	

\*Not anodized.



LOCKING SCREWS - Ø3.5 mm*			
Ref.	Description	Qty	
SOT3.5L10-ST	Locking screw - Ø3.5 mm - L 10 mm - STERILE	2	
SOT3.5L12-ST	Locking screw - Ø3.5 mm - L 12 mm - STERILE	1	
SOT3.5L14-ST	Locking screw - Ø3.5 mm - L 14 mm - STERILE	1	
SOT3.5L16-ST	Locking screw - Ø3.5 mm - L 16 mm - STERILE	1	
SOT3.5L18-ST	Locking screw - Ø3.5 mm - L 18 mm - STERILE	2	
SOT3.5L20-ST	Locking screw - Ø3.5 mm - L 20 mm - STERILE	2	
SOT3.5L22-ST	Locking screw - Ø3.5 mm - L 22 mm - STERILE	2	
SOT3.5L24-ST	Locking screw - Ø3.5 mm - L 24 mm - STERILE	2	

\*Blue anodized.



STANDARD CORTICAL SCREWS - Ø3.5 mm*			
Ref.	Description	Qty	
CT3.5L10-ST	Standard cortical screw - Ø3.5 mm - L 10 mm - STERILE	2	
CT3.5L12-ST	Standard cortical screw - Ø3.5 mm - L 12 mm - STERILE	1	
CT3.5L14-ST	Standard cortical screw - Ø3.5 mm - L 14 mm - STERILE	1	
CT3.5L16-ST	Standard cortical screw - Ø3.5 mm - L 16 mm - STERILE	1	
CT3.5L18-ST	Standard cortical screw - Ø3.5 mm - L 18 mm - STERILE	1	
CT3.5L20-ST	Standard cortical screw - Ø3.5 mm - L 20 mm - STERILE	2	
CT3.5L22-ST	Standard cortical screw - Ø3.5 mm - L 22 mm - STERILE	2	
CT3.5L24-ST	Standard cortical screw - Ø3.5 mm - L 24 mm - STERILE	2	

\* Light blue anodized.



### Also available on demand

Longer standard cortical screws (30 to 38 mm long) are available on demand.

To order, use the code CT3.5Lxx-ST and replace "xx" by the desired length. Example: "CT3.5L30-ST"

## > Removal and rescue kits

Sterile instruments

REMOVAL KITS		
Ref.	Description	Content
KIT-REMOVE-1	Removal kit for hexagonal stamp 2.0 mm	• 1x 2.0 mm hexagonal prehensor screwdriver
KIT-REMOVE-A	Removal kit for hexagonal stamp 2.5 mm	• 1x 2 in 1: 2.5 mm hexagonal prehensor screwdriver - Ø3.5 mm countersink

RESCUE KITS		
Ref.	Description	Content
KIT-RESCUE-2*	Rescue kit for Ø2.8 mm screws - Initial C & Initial A	• 1x Ø2.0 mm quick coupling drill bit - L 125 mm • 1x Ø2.0 mm threaded guide gauge for Ø2.8 mm screws
KIT-RESCUE-7	Rescue kit for Ø3.5 mm screws - Initial A	• 1x Ø2.7 mm quick coupling drill bit - L 125 mm • 1x Ø2.7 mm threaded guide gauge for Ø3.5 mm screws • 1x Handle for guide gauge • 1x Length gauge for Ø2.8 and Ø3.5 mm screws - L 10-60 mm • 2x Pins - Ø1.4 L120 mm



\* WARNING: When using a drill or locking guide gauge from a rescue kit, only use it in combination with the instruments of the same rescue kit to guarantee an accurate measure.



Supplemental instrumentation kits

## > Templates

Sterile templates

INITIAL A™ TEMPLATES			
Ref.	Description	Qty	
ANC607	Standard lateral fibula template plate - Sizes 1/2/3 - Right side and left side	1	
			
ANC659	Narrow lateral fibula template plate - Sizes 1/2	1	
			

This information is intended to demonstrate the Newclip Technics portfolio of medical devices. Always refer to the package insert, product label and/or user instructions before using any Newclip Technics product. These products must be handled and/or implanted by trained and qualified staff who have read the instructions before use. A surgeon must always rely on her or his own professional clinical judgement when deciding whether to use a particular product when treating a particular patient. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Newclip Technics representative if you have questions about the availability of Newclip Technics products in your area.

# KIT-AL2D




SINGLE USE KIT

STERILE R






MADE IN FRANCE


Example of kit content:



Implants material: Titanium TA6V - ISO 5832-3 / ASTM F136  
Degree of accuracy for devices with a measuring function: ± 0.8 mm



## Left Distal Fibula Lateral Standard Size 2



**NEWCLIP TECHNICS (HQ)**  
PA de la Lande Saint Martin  
45 rue des Garotières  
44115 Haute-Goulaine, France  
+33 (0)2 28 21 23 25  
orders@newcliptechnics.com  
www.newcliptechnics.com

**NEWCLIP TECHNICS IBERIA**  
Calle Frederic Mompou, 4b  
Sant Just Desvern  
08960 Barcelona, España  
+34 938 299 526  
contact@newclipiberia.com  
www.newcliptechnics.com

**NEWCLIP TECHNICS USA**  
Newclip USA  
642 Larkfield Center  
Santa Rosa CA 95403, USA  
+1 707 230 5078  
customerservice@newclipusa.com  
www.newclipusa.com

**NEWCLIP TECHNICS GERMANY**  
Newclip GmbH  
Pröllstraße 11, D-86157 Augsburg,  
Deutschland  
+49 (0)821 650 749 40  
info@newclipgmbh.com  
www.newclipgmbh.de

**NEWCLIP TECHNICS JAPAN**  
Newclip Technics Japan K.K.  
KKK Bldg. 502, 3-18-1 Asakusabashi  
Taito-Ku, Tokyo, 111-0053, Japan  
+81 (0)3 58 25 49 81  
www.newcliptechnics.com

**NEWCLIP TECHNICS AUSTRALIA**  
Newclip Australia  
3B/11 Donkin Street  
West End 4101, Australia  
+61 (0)2 81 886 110  
solutions@newclipaustralia.com  
www.newcliptechnics.com

Manufacturer: Newclip Technics - Brochure EN - Initial A - ED 07 - 08/21 - Medical device EC, class IIb - CE f039 SGS BE  
Read labeling and instructions before the use of Newclip Technics medical devices. Read labeling and instructions before the use of Newclip Technics medical devices.  
These products must be handled and/or implanted by trained and qualified staff who have read the instructions before use.

Non contractual pictures.