




NEWCLIP  
TECHNICS

---

## SURGICAL TECHNIQUE



---

### XPERT WRIST 2.4

DISTAL RADIUS VOLAR RIM PLATES


# DISTAL RADIUS VOLAR RIM PLATES

## SURGICAL TECHNIQUE - DISTAL RADIUS VOLAR RIM PLATES (PAGE 1/2)

**Example: volar rim plate (DETDVN1) for distal radius - Narrow head**




1. Position the plate onto the watershed line using the lateral lip of the plate.



2. Drill (ANC696) using the threaded guide gauge (ANC694) or the non-threaded bent guide gauge (ANC695) into the oblong hole.


Determine the screw length directly on the guide gauge or use the length gauge (ANC102).



3. Insert the Ø2.4 mm standard cortical screw (CT2.4Lxx) using the screwdriver (ANC575).

**NB:** in case of poor bone quality, insert a Ø2.4 mm locking screw (SDT2.4Lxx).

### OPTIONAL STEP



4. To ensure that the screws do not go into the joint, insert the pin (33.0214.120) into the radioulnar pin hole and verify its position by x-ray.

If necessary, remove the pin and readjust the plate positioning using the oblong hole.



5. Lock the threaded guide gauge (ANC694) in the radioulnar locking hole.

Determine the screw length directly on the guide gauge (ANC694) or use the length gauge (ANC102).

**NB:** if necessary you can modify the angulation using the polyaxial drill guide (ANC687) and the drill (ANC696). Then measure the depth using the length gauge (ANC102).




6. Insert a Ø2.4 mm locking screw (SDT2.4Lxx) using the screwdriver (ANC575).

# DISTAL RADIUS VOLAR RIM PLATES

## SURGICAL TECHNIQUE - DISTAL RADIUS VOLAR RIM PLATES (PAGE 2/2)

Example: evolar rim plate (DETDVN1) for distal radius - Narrow head



7. Repeat the last 2 steps for the remaining locking screws (SDT2.4Lxx).

### OPTIONAL STEP



8. In order to support the distal part, proceed in the same way as steps 5 and 6 (see page 14) for the monoaxial hole in the window.



### FINAL RESULT

#### Post-operative consideration

The plate positioning onto the watershed line may increase the risk of tendon injury. The surgeon should take this into consideration during subsequent follow-up of the patient. Plate removal post-healing is mandatory.

# IMPLANT REFERENCES

## DISTAL RADIUS VOLAR RIM PLATES

### DISTAL RADIUS VOLAR RIM PLATES

| Ref.      | Description                                                           |
|-----------|-----------------------------------------------------------------------|
| DETVN1 •  | Extra distal plate for distal radius - Narrow head - Size 1 - Left    |
| DETDVN1 • | Extra distal plate for distal radius - Narrow head - Size 1 - Right   |
| DETVS1 •  | Extra distal plate for distal radius - Standard head - Size 1 - Left  |
| DETDVS1 • | Extra distal plate for distal radius - Standard head - Size 1 - Right |
| DETVW1 •  | Extra distal plate for distal radius - Wide head - Size 1 - Left      |
| DETDW1 •  | Extra distal plate for distal radius - Wide head - Size 1 - Right     |

- Available in an Initial R Xpert single use kit (see below for more information)

### Initial R™ XPERT 2.4 single use kit

Newclip Technics also offers a single use sterile solution to treat hand and forearm fractures, osteotomies and arthrodeses: **Initial R™ Xpert 2.4**.

These kits are a range of single use kits with instrumentation and implants ready to use.

For more information, please refer to the Initial R Xpert brochure.

Please contact your NEWCLIP TECHNICS representative if you have questions about the availability of NEWCLIP TECHNICS products in your area.




# SCREW REFERENCES

## → Ø2.4 MM SCREWS

### Remark:


All implants are also available in a sterile version. An «-ST» is added to the end of the reference. Ex: «SDT2.4L08-ST»



### Ø2.4 MM LOCKING SCREWS\*

| Ref.      | Description                                        |
|-----------|----------------------------------------------------|
| SDT2.4L08 | Locking screw with conical head - Ø2.4 mm - L08 mm |
| SDT2.4L10 | Locking screw with conical head - Ø2.4 mm - L10 mm |
| SDT2.4L12 | Locking screw with conical head - Ø2.4 mm - L12 mm |
| SDT2.4L14 | Locking screw with conical head - Ø2.4 mm - L14 mm |
| SDT2.4L16 | Locking screw with conical head - Ø2.4 mm - L16 mm |
| SDT2.4L18 | Locking screw with conical head - Ø2.4 mm - L18 mm |
| SDT2.4L20 | Locking screw with conical head - Ø2.4 mm - L20 mm |
| SDT2.4L22 | Locking screw with conical head - Ø2.4 mm - L22 mm |
| SDT2.4L24 | Locking screw with conical head - Ø2.4 mm - L24 mm |
| SDT2.4L26 | Locking screw with conical head - Ø2.4 mm - L26 mm |
| SDT2.4L28 | Locking screw with conical head - Ø2.4 mm - L28 mm |
| SDT2.4L30 | Locking screw with conical head - Ø2.4 mm - L30 mm |

\* Non anodized



### Ø2.4 MM CORTICAL SCREWS\*

| Ref.     | Description                                |
|----------|--------------------------------------------|
| CT2.4L08 | Standard cortical screw - Ø2.4 mm - L08 mm |
| CT2.4L10 | Standard cortical screw - Ø2.4 mm - L10 mm |
| CT2.4L12 | Standard cortical screw - Ø2.4 mm - L12 mm |
| CT2.4L14 | Standard cortical screw - Ø2.4 mm - L14 mm |
| CT2.4L16 | Standard cortical screw - Ø2.4 mm - L16 mm |
| CT2.4L18 | Standard cortical screw - Ø2.4 mm - L18 mm |
| CT2.4L20 | Standard cortical screw - Ø2.4 mm - L20 mm |
| CT2.4L22 | Standard cortical screw - Ø2.4 mm - L22 mm |
| CT2.4L24 | Standard cortical screw - Ø2.4 mm - L24 mm |
| CT2.4L26 | Standard cortical screw - Ø2.4 mm - L26 mm |
| CT2.4L28 | Standard cortical screw - Ø2.4 mm - L28 mm |
| CT2.4L30 | Standard cortical screw - Ø2.4 mm - L30 mm |

\* Pink anodized

# INSTRUMENT REFERENCES

| INSTRUMENTS        |                                                                     |     |
|--------------------|---------------------------------------------------------------------|-----|
| Ref.               | Description                                                         | Qty |
| ANC102             | Length gauge                                                        | 1   |
| ANC350             | Ø4.5 mm AO quick coupling handle - Size 1                           | 1   |
| ANC503             | 150 mm reduction forceps                                            | 1   |
| ANC504             | 150 mm pointed reduction forceps                                    | 1   |
| ANC575             | T8 quick coupling screwdriver                                       | 2   |
| ANC578*            | Bending pliers                                                      | 2   |
| ANC687             | Polyaxial drill guide - SDT2.4 hole                                 | 2   |
| ANC694             | Ø1.8 mm threaded guide gauge for Ø2.4 mm screws                     | 2   |
| ANC695             | Ø1.8 mm non threaded bent guide gauge for Ø2.4 mm screws            | 1   |
| ANC696             | Ø1.8 mm quick coupling drill bit - L125 mm                          | 2   |
| ANC904** •         | MIS distal guide for distal radius - Narrow head - Left             | 1   |
| ANC905** •         | MIS distal guide for distal radius - Narrow head - Right            | 1   |
| ANC906** •         | MIS distal guide for distal radius - Standard head - Left           | 1   |
| ANC907** •         | MIS distal guide for distal radius - Standard head - Right          | 1   |
| ANC1061 ** •       | MIS distal guide for distal radius - Extra short - Right            | 1   |
| ANC1062** •        | MIS distal guide for distal radius - Extra short - Left             | 1   |
| ANC908** •         | Ø1.8 mm non threaded guide gauge                                    | 1   |
| ANC909** •         | Ø1.8 mm threaded guide gauge - MIS Xpert                            | 2   |
| ANC910** •         | T8 screwdriver with AO quick coupling system                        | 1   |
| ANC978             | Pin support for Ø1.4 mm pin - Long                                  | 1   |
| ANC979** •         | Pin support for Ø1.2 mm pin                                         | 1   |
| ANC1011            | Rim hook prehensor                                                  | 1   |
| ANC1023** •        | Ø1.2 mm short pin guide for distal radius arthroscopy               | 6   |
| ANC1024** •        | Ø1.2 mm long pin guide for distal radius arthroscopy                | 4   |
| ANC1025** ** •     | 2 in 1 instrument : T8 screwdriver Ø1.3 mm cannulated - Guide gauge | 1   |
| ANC1026** •        | Ø2.45 mm quick coupling drill bit - cannula Ø1.3 mm                 | 1   |
| ANC1136** •        | Pins guides support                                                 | 1   |
| 33.0212.120-MAR1 • | Pin Ø1.2 L120 mm                                                    | 8   |
| 33.0214.120        | Pin Ø1.4 L120 mm                                                    | 6   |
| 33.0218.080 •      | Pin Ø1.8 L80 mm                                                     | 2   |
| TD-111401-1.0NM-B  | Ø4.5 mm AO quick coupling handle with torque driver 1Nm             | 1   |

\*\* Available on demande

• Minimally Invasive Surgery

• Arthroscopic dedicated instruments

## REMOVAL KIT

If you have to remove XPERT WRIST 2.4 implants, make sure to order the **Newclip Technics** removal set, which includes the following instruments:

- ANC575: T8 quick coupling screwdriver
- ANC350: Ø4.5 mm AO quick coupling handle - Size 1
- KIT-REMOVE-2 : sterile T8 prehensor screwdriver for Initial R Xpert kit.

## \*BENDING PLIERS


The bending pliers must not be used with the XPERT WRIST 2.4 - VOLAR PLATE RANGE

**Indications :** The implants of the Xpert Wrist 2.4 range are intended for hand and forearm fractures, osteotomies and arthrodeses in adults.

### Contraindications:

- Serious vascular deterioration, bone devitalization.
- Pregnancy.
- Acute or chronic local or systemic infections.
- Lack of musculo-cutaneous cover, severe vascular deficiency affecting the concerned area.
- Insufficient bone quality preventing a good fixation of the implants into the bone.
- Muscular deficit, neurological deficiency or behavioral disorders, which could submit the implant to abnormal mechanical strains.
- Allergy to one of the materials used or sensitivity to foreign bodies.
- Serious problems of non-compliance, mental or neurological disorders, failure to follow post-operative care recommendations.
- Unstable physical and/or mental condition.

# KIT DESCRIPTION



The information presented in this brochure is intended to demonstrate a NEWCLIP TECHNICS product. Always refer to the package insert, product label and/or user instructions before using any NEWCLIP TECHNICS product. Surgeons must always rely on their own clinical judgment when deciding which products and techniques to use with their patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your NEWCLIP TECHNICS representative if you have questions about the availability of NEWCLIP TECHNICS products in your area.



# NEWCLIP TECHNICS

INNOVATION MEANS MOTION

#### NEWCLIP TECHNICS

PA de la Lande Saint Martin  
45 rue des Garottières  
44115 Haute Goulaine, France  
+33 (0)2 28 21 23 25  
orders@newcliptechnics.com  
www.newcliptechnics.com

#### NEWCLIP TECHNICS USA

Newclip USA  
340 Tesconi Circle Suite A  
Santa Rosa CA 95403, USA  
+1 707 230 5078  
customerservice@newclipusa.com  
www.newclipusa.com

#### NEWCLIP TECHNICS GERMANY

Newclip GmbH  
Pröllstraße 11  
D-86157 Augsburg, Deutschland  
+49 (0)821 650 749 40  
info@newclipgmbh.com  
www.newclipgmbh.de

#### NEWCLIP TECHNICS JAPAN

Newclip Technics Japan K.K.  
KKK Bldg. 502, 3-18-1 Asakusabashi  
Taito-Ku, Tokyo, 111-0053, Japan  
+81 (0)3 58 25 49 81  
Fax: +81 (0)3 58 25 49 86  
www.newcliptechnics.com

#### NEWCLIP TECHNICS AUSTRALIA

Newclip Australia  
3B/11 Donkin Street  
West End 4101, Australia  
+61 (0)2 81 886 110  
solutions@newclipaustralia.com  
www.newcliptechnics.com

MADE IN FRANCE