

Initial A™: Single use distal fibula sterile kits

Initial A™ - Fibula & Syndesmosis kits content

Intended purpose :


The implants of the Initial A™ range are intended for the fixation of fractures, osteotomies and pseudarthroses of the distal and the diaphyseal fibula, the distal tibia and for the syndesmosis repair in adults.

Contraindications:


- Pregnancy.
- Acute or chronic local or systemic infections.
- Allergy to one of the materials used or sensitivity to foreign bodies.



> Initial A™ - Fibula kit




> Initial A™ - Syndesmosis kit



NEWCLIP
TECHNICS




Initial A™ - Fibula kits references: Implants



Precontoured lateral plates
with polyaxial locking fixation
Standard - Size 1

KIT-AL1G or KIT-AL1D* Distal fibula kit - Standard - Left or right - Size 1		
Ref.	Description	Qty
RTxLS1	Lateral plate - Distal fibula - Standard - Size 1 (left or right)	1
SDT2.8L10	Ø2.8 mm locking screw - L10 mm	1
SDT2.8L12	Ø2.8 mm locking screw - L12 mm	1
SDT2.8L14	Ø2.8 mm locking screw - L14 mm	2
SDT2.8L16	Ø2.8 mm locking screw - L16 mm	2
SDT2.8L18	Ø2.8 mm locking screw - L18 mm	2
SOT3.5L12	Ø3.5 mm locking screw - L12 mm	1
SOT3.5L14	Ø3.5 mm locking screw - L14 mm	1
SOT3.5L16	Ø3.5 mm locking screw - L16 mm	1
SOT3.5L18	Ø3.5 mm locking screw - L18 mm	1
CT3.5L14	Ø3.5 mm non-locking screw - L14 mm	1
CT3.5L16	Ø3.5 mm non-locking screw - L16 mm	1
CT3.5L20	Ø3.5 mm non-locking screw - L20mm	1
CT3.5L24	Ø3.5 mm non-locking screw - L24 mm	1



Precontoured lateral plates
with polyaxial locking fixation
Standard - Size 2

KIT-AL2G or KIT-AL2D* Distal fibula kit - Standard - Left or right - Size 2		
Ref.	Description	Qty
RTxLS2	Lateral plate for distal fibula - Standard - Size 2 (left or right)	1
SDT2.8L10	Ø2.8 mm locking screw - L10 mm	1
SDT2.8L12	Ø2.8 mm locking screw - L12 mm	1
SDT2.8L14	Ø2.8 mm locking screw - L14 mm	2
SDT2.8L16	Ø2.8 mm locking screw - L16 mm	2
SDT2.8L18	Ø2.8 mm locking screw - L18 mm	2
SOT3.5L12	Ø3.5 mm locking screw - L12 mm	2
SOT3.5L14	Ø3.5 mm locking screw - L14 mm	2
SOT3.5L16	Ø3.5 mm locking screw - L16 mm	1
SOT3.5L18	Ø3.5 mm locking screw - L18 mm	1
CT3.5L12	Ø3.5 mm non-locking screw - L12 mm	1
CT3.5L14	Ø3.5 mm non-locking screw - L14 mm	1
CT3.5L16	Ø3.5 mm non-locking screw - L16mm	1
CT3.5L18	Ø3.5 mm non-locking screw - L18 mm	1
CT3.5L22	Ø3.5 mm non-locking screw - L22 mm	1




Precontoured lateral plates
with polyaxial locking fixation
Standard - Size 3

KIT-AL3G or KIT-AL3D* Distal fibula kit - Standard - Left or right - Size 3		
Ref.	Description	Qty
RTxLS3	Lateral plate - Distal fibula - Standard - Size 3 (left or right)	1
SDT2.8L10	Ø2.8 mm locking screw - L10 mm	1
SDT2.8L12	Ø2.8 mm locking screw - L12 mm	1
SDT2.8L14	Ø2.8 mm locking screw - L14 mm	2
SDT2.8L16	Ø2.8 mm locking screw - L16 mm	2
SDT2.8L18	Ø2.8 mm locking screw - L18 mm	2
SOT3.5L12	Ø3.5 mm locking screw - L12 mm	3
SOT3.5L14	Ø3.5 mm locking screw - L14 mm	2
SOT3.5L16	Ø3.5 mm locking screw - L16 mm	1
SOT3.5L18	Ø3.5 mm locking screw - L18 mm	1
CT3.5L12	Ø3.5 mm non-locking screw - L12 mm	1
CT3.5L14	Ø3.5 mm non-locking screw - L14 mm	2
CT3.5L16	Ø3.5 mm non-locking screw - L16mm	1
CT3.5L18	Ø3.5 mm non-locking screw - L18 mm	1
CT3.5L22	Ø3.5 mm non-locking screw - L22 mm	1



Precontoured lateral plate
with polyaxial locking fixation
Narrow - Size 1

KIT-AL1S Distal fibula kit - Narrow symmetrical - Size 1		
Ref.	Description	Qty
RTSLN1	Lateral plate - Distal fibula - Narrow - Symmetrical - Size 1	1
SDT2.8L14	Ø2.8 mm locking screw - L14 mm	1
SDT2.8L16	Ø2.8 mm locking screw - L16 mm	2
SDT2.8L18	Ø2.8 mm locking screw - L18 mm	1
SOT3.5L12	Ø3.5 mm locking screw - L12 mm	2
SOT3.5L14	Ø3.5 mm locking screw - L14 mm	2
SOT3.5L16	Ø3.5 mm locking screw - L16 mm	1
CT3.5L14	Ø3.5 mm non-locking screw - L14 mm	1
CT3.5L16	Ø3.5 mm non-locking screw - L16 mm	1
CT3.5L22	Ø3.5 mm non-locking screw - L22 mm	1



Precontoured lateral plate
with polyaxial locking fixation
Narrow - Size 2

KIT-AL2S Distal fibula kit - Narrow symmetrical - Size 2		
Ref.	Description	Qty
RTSLN2	Lateral plate - Distal fibula - Narrow - Symmetrical - Size 2	1
SDT2.8L14	Ø2.8 mm locking screw - L14 mm	1
SDT2.8L16	Ø2.8 mm locking screw - L16 mm	2
SDT2.8L18	Ø2.8 mm locking screw - L18 mm	1
SOT3.5L12	Ø3.5 mm locking screw - L12 mm	2
SOT3.5L14	Ø3.5 mm locking screw - L14 mm	2
SOT3.5L16	Ø3.5 mm locking screw - L16 mm	2
CT3.5L12	Ø3.5 mm non-locking screw - L12 mm	1
CT3.5L14	Ø3.5 mm non-locking screw - L14 mm	1
CT3.5L16	Ø3.5 mm non-locking screw - L16 mm	1
CT3.5L22	Ø3.5 mm non-locking screw - L22 mm	1

* Each kit includes one plate (left or right)
NB : Additional sterile screws can be ordered separately .



NEWCLIP TECHNICS (HQ)
45 rue des Garottières
44115 Haute Goulaine, France
+33 (0)2 28 21 23 25
orders@newcliptechnics.com
www.newcliptechnics.com

NEWCLIP TECHNICS USA
340 Tesconi Circle Suite A
Santa Rosa CA 95401, USA
+1 707 230 5078
customerservice@newclipusa.com
www.newcliptechnics.com

NEWCLIP TECHNICS GERMANY
Pröllstraße 11
D-86157 Augsburg, Deutschland
+49 (0)821 650 749 40
info@newclipgmbh.com
www.newcliptechnics.com

NEWCLIP TECHNICS AUSTRALIA
38/11 Donkin Street
West End 4101, Australia
+61 (0)2 81 886 110
solutions@newclipaustralia.com
www.newcliptechnics.com

NEWCLIP TECHNICS JAPAN
KKK Bldg 502, 3-18-1 Asakusabashi
Taisho-Ku, Tokyo, 111-0053, Japan
+81 (0)3 58 25 49 81
Fax: +81 (0)3 58 25 49 86
www.newcliptechnics.com

NEWCLIP TECHNICS BELGIUM
Av. du Bourg Etienne
Demunter 5
1090 Jette, Belgium
+32 (0)9 23 11 33 1
www.newcliptechnics.fr

NEWCLIP TECHNICS IBERIA
Calle Frederic Mompou, 4b
Sant Just Desvern, 08960 Barcelona, España
+34 938 299 526
contact@newclipiberia.com
www.newcliptechnics.com

Initial A™ - Syndesmosis Kit references

INITIAL A™ - SYNDESMOSIS KIT

Ref.	Description
KIT-AMS	Syndesmosis kit

ASSOCIATED STERILE SCREWS - AVAILABLE SEPARATELY

Ref.	Description	Qty
CT3.5L40-ST	Ø3.5 mm non-locking screw - L40 mm - STERILE	2
CT3.5L45-ST	Ø3.5 mm non-locking screw - L45 mm - STERILE	3
CT3.5L50-ST	Ø3.5 mm non-locking screw - L50 mm - STERILE	3
CT3.5L55-ST	Ø3.5 mm non-locking screw - L55 mm - STERILE	3
CT3.5L60-ST	Ø3.5 mm non-locking screw - L60 mm - STERILE	2
CT3.5L65-ST	Ø3.5 mm non-locking screw - L65 mm - STERILE	2
QT4.0L40-ST	Ø4.0 mm lag screw - L40 mm - STERILE	3
QT4.0L45-ST	Ø4.0 mm lag screw - L45 mm - STERILE	3
QT4.0L50-ST	Ø4.0 mm lag screw - L50 mm - STERILE	3
QT4.0L55-ST	Ø4.0 mm lag screw - L55 mm - STERILE	3
QT4.0L60-ST	Ø4.0 mm lag screw - L60 mm - STERILE	3
QT4.0L65-ST	Ø4.0 mm lag screw - L65 mm - STERILE	3