

# Initial C™: Single use clavicle sterile kits

## Initial C™ - Clavicle kits contents




### Intended purpose:

The implants of the Initial C™ range are dedicated to the fixation of fractures, mal-union, non-unions, and osteotomies of the clavicle in adults.

### Contraindications:

- Pregnancy.
- Acute or chronic local or systemic infections.
- Allergy to one of the materials used or sensitivity to foreign bodies.



NEWCLIP-TECHNICS



# Initial C™: Single use clavicle sterile kits

## Initial C™ - Clavicle kits contents

### > Initial C™ kits references - Midshaft and lateral midshaft plates

#### MIDSHAFT AND MIDSHAFT- LATERAL PLATES

Ref.	Description
KIT-CM2D	Clavicle kit - Superior - Midshaft - Size 2 - Right
KIT-CM2G	Clavicle kit - Superior - Midshaft - Size 2 - Left
KIT- CML1D	Clavicle kit - Superior - Lateral midshaft - Size 1 - Right
KIT- CML1G	Clavicle kit - Superior - Lateral midshaft - Size 1 - Left
KIT- CML2D	Clavicle kit - Superior - Lateral midshaft - Size 2 - Right
KIT- CML2G	Clavicle kit - Superior - Lateral midshaft - Size 2 - Left

#### BENDABLE MIDSHAFT PLATES

Ref.	Description
KIT-CBM1D	Clavicle kit - Superior - Midshaft bendable - Size 1 - Right
KIT-CBM1G	Clavicle kit - Superior - Midshaft bendable - Size 1 - Left
KIT-CBM2D	Clavicle kit - Superior - Midshaft bendable - Size 2 - Right
KIT-CBM2G	Clavicle kit - Superior - Midshaft bendable - Size2 - Left

#### INITIAL C™ KITS CONTENT - IMPLANTS

		QUANTITY PER KITS					
	Ref.	Description	KIT-CBM1D/G	KIT-CM2D/G	KIT-CBM2D/G	KIT-CML1D/G	KIT-CML2D/G
PLATES	CBTDM1D or CBTGM1D	Midshaft bendable clavicle plate - Superior - Size 1 - Right or Left	1	-	-	-	-
	CTDM2D or CTGM2D	Midshaft clavicle plate - Superior - Size 2 - Right or Left	-	1	-	-	-
	CBTDM2D or CBTGM2D	Midshaft bendable clavicle plate - Superior - Size 2 - Right or Left	-	-	1	-	-
	CTDML1D or CTGML1D	Lateral midshaft clavicle plate - Superior - Size 1 - Right or Left	-	-	-	1	-
	CTDML2D or CTGML2D	Lateral midshaft clavicle plate - Superior - Size 2 - Right or Left	-	-	-	-	1
NON-LOCKING SCREWS Ø3.5 MM	CT3.5L12D	Ø3.5 mm non-locking screw - L12 mm	1	-	1	1	-
	CT3.5L14D	Ø3.5 mm non-locking screw - L14 mm	1	2	1	1	2
	CT3.5L16D	Ø3.5 mm non-locking screw - L16 mm	1	2	1	1	2
	CT3.5L18D	Ø3.5 mm non-locking screw - L18 mm	1	-	1	1	-
LOCKING SCREWS Ø3.5 MM	SOT3.5L12D	Ø3.5 mm locking screw - L12 mm	1	1	2	1	2
	SOT3.5L14D	Ø3.5 mm locking screw - L14 mm	2	3	4	2	4
	SOT3.5L16D	Ø3.5 mm locking screw - L16 mm	2	3	3	2	3
	SOT3.5L18D	Ø3.5 mm locking screw - L18 mm	1	1	2	1	-

# Initial C™: Single use clavicle sterile kits

## Initial C™ - Clavicle kits contents

INITIAL C™ KITS CONTENT - IMPLANTS			QUANTITY PER KITS		
	Ref.	Description	KIT-CL1D/G	KIT-CL2D/G	KIT-CBL2D/G
PLATES	CTDL1D	Lateral clavicle plate - Superior - Size 1 - Right or Left	1	-	-
	CTGL1D				
	CTDL2D	Lateral clavicle plate - Superior - Size 2 - Right or Left	-	1	-
	CTGL2D				
CBTDL2D	Lateral bendable clavicle plate - Superior - Size 2 - Right or Left	-	-	1	
CBTGL2D					
LOCKING SCREWS Ø2.8 MM	SDT2.8L10D	Ø2.8 mm locking screw - L10 mm	1	1	1
	SDT2.8L12D	Ø2.8 mm locking screw - L12 mm	2	2	2
	SDT2.8L14D	Ø2.8 mm locking screw - L14 mm	2	2	2
	SDT2.8L16D	Ø2.8 mm locking screw - L16 mm	2	2	2
	SDT2.8L18D	Ø2.8 mm locking screw - L18 mm	1	1	1
NON-LOCKING SCREWS Ø3.5 MM	CT3.5L12D	Ø3.5 mm non-locking screw - L12 mm	1	1	1
	CT3.5L14D	Ø3.5 mm non-locking screw - L14 mm	1	1	1
	CT3.5L16D	Ø3.5 mm non-locking screw - L16 mm	1	1	1
	CT3.5L18D	Ø3.5 mm non-locking screw - L18 mm	0	1	1
LOCKING SCREWS Ø3.5 MM	SOT3.5L12D	Ø3.5 mm locking screw - L12 mm	1	1	1
	SOT3.5L14D	Ø3.5 mm locking screw - L14 mm	2	2	2
	SOT3.5L16D	Ø3.5 mm locking screw - L16 mm	1	2	2
	SOT3.5L18D	Ø3.5 mm locking screw - L18 mm	1	1	3

## > Initial C™ kits references - Lateral plates

LATERAL PLATES	
Ref.	Description
KIT-CL1D	Clavicle kit - Superior - Lateral - Size 1 - Right
KIT-CL1G	Clavicle kit - Superior - Lateral - Size 1 - Left
KIT-CL2D	Clavicle kit - Superior - Lateral - Size 2 - Right
KIT-CL2G	Clavicle kit - Superior - Lateral - Size 2 - Left
KIT-CBL2D	Clavicle kit - Superior - Lateral bendable - Size 2 - Right
KIT-CBL2G	Clavicle kit - Superior - Lateral bendable - Size 2 - Left

This information is intended to demonstrate the Newclip Technics portfolio of medical devices. Always refer to the package insert, product label and/or user instructions before using any Newclip Technics product. These products must be handled and/or implanted by trained and qualified staff who have read the instructions before use. A surgeon must always rely on her or his own professional clinical judgement when deciding whether to use a particular product when treating a particular patient. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Newclip Technics representative if you have questions about the availability of Newclip Technics products in your area.



### NEWCLIP TECHNICS (HQ)

45 rue des Garottières  
44115 Haute Goulaine, France  
+33 (0)2 28 21 23 25  
orders@newcliptechnics.com  
www.newcliptechnics.com

### NEWCLIP TECHNICS USA

Newclip USA LLC  
340 Tesconi Circle, Suite A  
Santa Rosa CA 95401 USA  
+1 707 230 5078  
customerservice@newclipusa.com  
www.newcliptechnics.com

### NEWCLIP TECHNICS IBERIA

Calle Frederic Mompou 4b  
Sant Just Desvern  
08960 Barcelona, España  
+34 938 299 526  
contact@newclipiberia.com  
www.newcliptechnics.com

### NEWCLIP TECHNICS GERMANY

Newclip GmbH  
Pröllstraße 11  
D-86157 Augsburg, Deutschland  
+49 (0)821 650 749 40  
info@newclipgmbh.com  
www.newcliptechnics.com

### NEWCLIP TECHNICS AUSTRALIA

Newclip Australia  
3B/11 Donkin Street  
West End 4101, Australia  
+61 (0)2 81 886 110  
solutions@newclipaustralia.com  
www.newcliptechnics.com

### NEWCLIP TECHNICS BELGIUM

Newclip Belgium  
Detbystraat 43  
9051 Sint-Denijs-Westrem  
+32 (0)9 2311331  
contact@newclipbelgium.com  
www.newcliptechnics.com