

Initial F™-Lapidus: single use sterile kits for plantar lapidus arthrodesis

Kit content


Example : Kit LV1D - CQ4.0

> Intended purpose

The implants of the Initial F™ range are intended for arthrodeses, fractures and osteotomies fixation and revision surgeries of the foot in adults.

> Contraindications

- Pregnancy
- Acute or chronic, local or systemic infections.
- Allergy to one of the materials used or sensitivity to foreign bodies.




Initial F™ Lapidus: single use sterile kits for plantar lapidus arthrodesis

Initial F™ Lapidus - Kits references


> Kits content - Implants



Plantar Lapidus arthrodesis plate - Narrow - Size 0 - Left or Right



Plantar Lapidus arthrodesis plate - Narrow - Size 1 - Left or Right



Plantar Lapidus arthrodesis plate - Narrow - Size 2 - Left or Right

KIT-LV0G-CQ4.0 or KIT-LV0D-CQ4.0*			
Foot kit - Lapidus Arthrodesis - Plantar - Size 0 - Left + for Ø4.0 mm screws			
Ref.	Description	Length	Qty
FLTvx0	Plantar Lapidus arthrodesis plate - Narrow - Size 0 - Left or Right	37.0	1
RLT3.5L20	Ø3.5 mm non-locking screw - L20 mm		1
RLT3.5L26	Ø3.5 mm non-locking screw - L26 mm		1
RLT3.5L28	Ø3.5 mm non-locking screw - L28mm		1
SLT3.5L14	Ø3.5 mm locking screw - L14 mm		1
SLT3.5L16	Ø3.5 mm locking screw - L16 mm		1
SLT3.5L20	Ø3.5 mm locking screw - L20 mm		1
SLT3.5L26	Ø3.5 mm locking screw - L26 mm		1

KIT-LV1G-CQ4.0 or KIT-LV1D-CQ4.0*			
Foot kit - Lapidus Arthrodesis - Plantar - Size 1 - Left + for Ø4.0 mm screws			
Ref.	Description	Length	Qty
FLTvx1	Plantar Lapidus arthrodesis plate - Narrow - Size 1 - Left or Right	41.8	1
RLT3.5L20	Ø3.5 mm non-locking screw - L20 mm		1
RLT3.5L26	Ø3.5 mm non-locking screw - L26 mm		1
RLT3.5L28	Ø3.5 mm non-locking screw - L28mm		1
SLT3.5L14	Ø3.5 mm locking screw - L14 mm		1
SLT3.5L16	Ø3.5 mm locking screw - L16 mm		1
SLT3.5L20	Ø3.5 mm locking screw - L20 mm		1
SLT3.5L26	Ø3.5 mm locking screw - L26 mm		1

KIT-LV2G-CQ4.0 or KIT-LV2D-CQ4.0*			
Foot kit - Lapidus Arthrodesis - Plantar - Size 2 - Left + for Ø4.0 mm screws			
Ref.	Description	Length	Qty
FLTvx2	Plantar Lapidus arthrodesis plate - Narrow - Size 2 - Left or Right	49.8	1
RLT3.5L20	Ø3.5 mm non-locking screw - L20 mm		1
RLT3.5L26	Ø3.5 mm non-locking screw - L26 mm		1
RLT3.5L28	Ø3.5 mm non-locking screw - L28mm		1
SLT3.5L14	Ø3.5 mm locking screw - L14 mm		2
SLT3.5L16	Ø3.5 mm locking screw - L16 mm		1
SLT3.5L20	Ø3.5 mm locking screw - L20 mm		1
SLT3.5L26	Ø3.5 mm locking screw - L26 mm		1

> Kits content - Instruments

STERILE INSTRUMENTATION FOR PLANTAR LAPIDUS PLATE	
Description	Qty
5.8 mm single use handle	1
Handle for guide gauge	1
Length gauge	1
Ø2.7 mm threaded guide gauge for Ø3.5 mm screws	1
Ø2.7 mm quick coupling drill bit - L 125 mm	1
T8 quick coupling screwdriver	1
Pin Ø1.6 - L100 mm	1
Pin Ø1.2 - L170 mm	2

STERILE INSTRUMENTATION FOR Ø4.0 MM CANNULATED SCREWS*	
Description	Qty
Length gauge for pin Ø1.3 mm	1
2 in 1: 2.5 mm hexagonal screwdriver - Ø6.0 mm countersink	1
Ø2.9 mm drill bit - cannula 1.4 mm - L 120 mm - AO Ø4.5 mm quick coupling	1
Pin Ø1.3 - L140 mm	3

* The Ø4.0 mm cannulated headed screws (Ref. H1.4QT4.0xxL-ST) are available in supplemental sterile screw caddy. The available range is from 26 mm to 48 mm with an incremental of 2 mm.



> Additional instrumentation kits

INITIAL F™ - LAPIDUS STERILE TEMPLATES	
Ref.	Description
ANC958	Template for foot kit - Lapidus Arthrodesis - Plantar - Sizes 0-1-2 - Left & Right (KIT-LVxx)

> Removal and Rescue kits

INITIAL F™ - REMOVAL AND RESCUE KITS	
Ref.	Description
KIT-RESCUE-6	Rescue kit for Initial F lapidus
Kit-REMOVE-A	Removal kit for hexagonal stamp 2.5 mm
KIT-REMOVE-2	Removal kit for T8 hexalobe



NEWCLIP TECHNICS
PA de la Lande Saint Martin - 45 rue des Garotières
44115 Haute Goulaine (France)
Tél. : +33 (0)2 28 21 23 25 - Fax : +33 (0)2 40 63 68 37
orders@newcliptechnics.com
www.newcliptechnics.com

NEWCLIP USA
642 Larkfield Center
Santa Rosa CA 95403, USA
Phone : +1 707 230 5078
customerservice@newclipusa.com
www.newcliptechnics.com

NEWCLIP TECHNICS GmbH
Wankelstraße 1B,
D-86391 Stadtbergen, Deutschland
Phone : +49 (0)821 650 749 40
info@newclipgmbh.com
www.newcliptechnics.com

NEWCLIP TECHNICS Australia
Forest Central Business Park
Building 10 Level 1, 49 Frenchs Forest
Rd, Frenchs Forest 2086, Australia
Phone: +61 (0)2 81 886 110
solutions@newclipaustralia.com
www.newcliptechnics.com

The information presented in this brochure is intended to demonstrate a Newclip Technics product. Always refer to the package insert, product label and/or user instructions before using any Newclip Technics product. Surgeons must always rely on their own clinical judgment when deciding which products and techniques to use with their patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Newclip Technics representative if you have questions about the availability of Newclip Technics products in your area.