

Initial F™ SiTa: Single use sterile kit for subtalar arthroereisis screws

Initial F™ SiTa- Kit content




> Intended purpose

The subtalar screws of the SiTa range are intended for hyperpronated foot surgery.

The Initial F™ SiTa kits consist of sets of instruments to be used only in orthopaedic surgery, when Newclip Technics subtalar screws are implanted or explanted.

> Contraindications

- Pregnancy.
- Acute or chronic local or systemic infections.
- Allergy to one of the materials used or sensitivity to foreign bodies.



NEWCLIP-TECHNICS

F
INITIAL SiTa

Initial F™ SiTa- References

Initial F™ SiTa kit

INITIAL F™ SITA KIT - INSTRUMENTATION CONTENT				
Ref.	Description	Content	Qty	Colour
KIT-SCY	Instrumentation kit for subtalar arthroereisis screws	Dilator for subtalar arthroereisis screw	1	-
		Screwdriver for subtalar arthroereisis screw	1	-
		Template for subtalar arthroereisis screw H2.1YT8.0L14D	1	●
		Template for subtalar arthroereisis screw H2.1YT9.0L14D	1	●
		Template for subtalar arthroereisis screw H2.1YT10.0L14D	1	●
		Template for subtalar arthroereisis screw H2.1YT11.0L16D	1	●
		Template for subtalar arthroereisis screw H2.1YT12.0L16D	1	●
		Pin Ø2.0 mm L250 mm	2	-

Implants

Subtalar arthroereisis screws available separately in a sterile version (S-box).

STERILE IMPLANTS		
Ref.	Description	Anodized screw colour
H2.1YT8.0L14D-ST	Ø8.0 mm subtalar arthroereisis screw - cannula Ø2.1 - L14 mm - STERILE	●
H2.1YT9.0L14D-ST	Ø9.0 mm subtalar arthroereisis screw - cannula Ø2.1 - L14 mm - STERILE	●
H2.1YT10.0L14D-ST	Ø10.0 mm subtalar arthroereisis screw - cannula Ø2.1 - L14 mm - STERILE	●
H2.1YT11.0L16D-ST	Ø11.0 mm subtalar arthroereisis screw - cannula Ø2.1 - L16 mm - STERILE	●
H2.1YT12.0L16D-ST	Ø12.0 mm subtalar arthroereisis screw - cannula Ø2.1 - L16 mm - STERILE	●

Removal kit

Sterile instruments available separately.

REMOVAL KIT			
Ref.	Description	Content	Qty
KIT-REMOVE-7	Removal kit for subtalar arthroereisis screw	- Screwdriver for subtalar arthroereisis screw	1
		- Pin Ø2.0 mm L250 mm	1



CORRESPONDING TABLE	
The Initial F™ SiTa kit includes 5 different size templates in a radiopaque plastic. They correspond to 5 associated subtalar arthroereisis screws. The template allows to select the appropriate size of implant for the desired correction.	
TEMPLATES	IMPLANTS
	 H2.1YT8.0L14D-ST
	 H2.1YT9.0L14D-ST
	 H2.1YT10.0L14D-ST
	 H2.1YT11.0L16D-ST
	 H2.1YT12.0L16D-ST

CAUTION: Do not use template to widen the sinus tarsi.

This information is intended to demonstrate the Newclip Technics portfolio of medical devices. Always refer to the package insert, product label and/or user instructions before using any Newclip Technics product. These products must be handled and/or implanted by trained and qualified staff who have read the instructions before use. A surgeon must always rely on her or his own professional clinical judgement when deciding whether to use a particular product when treating a particular patient. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Newclip Technics representative if you have questions about the availability of Newclip Technics products in your area.



NEWCLIP TECHNICS (HQ)

45 rue des Garotières
44115 Haute-Goulaine, France
+33 (0)2 28 21 23 25
orders@newcliptechnics.com
www.newcliptechnics.com

NEWCLIP TECHNICS IBERIA

Calle Frederic Mompou, 4b
Sant just Desvern
08960 Barcelona, España
+34 937 503 421
contact@newclipiberia.com
www.newcliptechnics.com

NEWCLIP TECHNICS USA

Newclip USA
340 Tesconi Circle, Suite A
Santa Rosa CA95401, USA
+1 707 230 5078
customerservice@newclipusa.com
www.newcliptechnics.com

NEWCLIP TECHNICS GERMANY

Newclip GmbH
Pröllstraße 11, D-86157 Augsburg,
Deutschland
+49 (0)821 650 749 40
info@newclipgmbh.com
www.newcliptechnics.com

NEWCLIP TECHNICS JAPAN

Newclip Technics Japan K.K.
KKK Bldg. 502, 3-18-1 Asakusabashi
Taito-Ku, Tokyo, 111-0053, Japan
+81 (0)3 58 25 49 81
www.newcliptechnics.com

NEWCLIP TECHNICS AUSTRALIA

Newclip Australia
3B/11 Donkin Street
West End 4101, Australia
+61 (0)2 81 886 110
solutions@newclipaustralia.com
www.newcliptechnics.com

NEWCLIP BELGIUM

Newclip Belgium
Av. du Bourg Etienne
Demunter 5
1090 Jette, Belgium
Phone +32 (0)9 2311331
contact@newclipbelgium.com
www.newcliptechnics.com