

Initial R Xpert™ 2.4 Fragment Specific: Single use distal radius and ulna sterile kits

Initial R Xpert 2.4 Fragment Specific - kits content



Example: KIT-INSTRUM-1 (on the left) and KIT-WHE2G (on the right)

> Intended purpose

The implants of the Initial R™ Xpert range are intended for fixation of hand and forearm fractures, osteotomies and arthrodeses in adults.


> Contraindications

- Pregnancy.
- Acute or chronic local or systemic infections.
- Allergy to one of the materials used or sensitivity to foreign bodies.



T8 prehensor screwdriver

Non threaded polyaxial drill guide



Length gauge

Ø1.8 mm threaded guide gauge


Ø1.8 mm quick coupling drill bit - L125 mm

Pins - Ø1.4 L120 mm (x4)



CT2.4Lxx
Ø2.4 mm non-locking screws
Pink anodized

SDT2.4Lxx
Ø2.4 mm locking screws
Non anodized



Fragment specific plate



NEWCLIP-TECHNICS

INITIAL
Xpert 2.4
Fragment
Specific

Initial R Xpert™ 2.4 Fragment Specific - Kit references



Dorso-medial kit – Distal radius – Extra Narrow – Size 2 – Right or Left

KIT-WRNS2x

KIT-WRNS2G OR KIT-WRNS2D* Dorso-medial kit - Implants only - Distal radius - Extra Narrow - Size 2 - Left or Right		
Ref.	Description	Qty
DTxRNS2	Dorso-medial plate - Distal radius - Extra Narrow - Size 2 - Left or Right	1
SDT2.4L12	Ø2.4 mm locking screw - L12 mm	2
SDT2.4L14	Ø2.4 mm locking screw - L14 mm	2
SDT2.4L16	Ø2.4 mm locking screw - L16 mm	2
SDT2.4L18	Ø2.4 mm locking screw - L18 mm	2
SDT2.4L20	Ø2.4 mm locking screw - L20 mm	2
SDT2.4L22	Ø2.4 mm locking screw - L22 mm	1
CT2.4L14	Ø2.4 mm non-locking screw - L14 mm	1
CT2.4L16	Ø2.4 mm non-locking screw - L16 mm	1
CT2.4L18	Ø2.4 mm non-locking screw - L18 mm	1



Latero-dorsal radial column kit – Distal radius – Size 2 – Right or Left

KIT-WEF2x

KIT-WEF2G OR KIT-WEF2D* Latero-dorsal radial column kit - Implants only - Distal radius - Size 2 - Left or Right		
Ref.	Description	Qty
DETxF2	Latero-dorsal radial column plate - Distal radius - Size 2 - Left or Right	1
SDT2.4L12	Ø2.4 mm locking screw - L12 mm	2
SDT2.4L14	Ø2.4 mm locking screw - L14 mm	2
SDT2.4L16	Ø2.4 mm locking screw - L16 mm	2
SDT2.4L18	Ø2.4 mm locking screw - L18 mm	2
SDT2.4L20	Ø2.4 mm locking screw - L20 mm	2
SDT2.4L22	Ø2.4 mm locking screw - L22 mm	1
CT2.4L14	Ø2.4 mm non-locking screw - L14 mm	1
CT2.4L16	Ø2.4 mm non-locking screw - L16 mm	1
CT2.4L18	Ø2.4 mm non-locking screw - L18 mm	1



Dorso-medial kit – Distal radius – Size 2 – Right or Left

KIT-WR2x


KIT-WR2G OR KIT-WR2D* Dorso-medial kit - Implants only - Distal radius - Size 2 - Left or Right		
Ref.	Description	Qty
DTxR2	Dorso-medial plate - Distal radius - Size 2 - Left or Right	1
SDT2.4L12	Ø2.4 mm locking screw - L12 mm	2
SDT2.4L14	Ø2.4 mm locking screw - L14 mm	2
SDT2.4L16	Ø2.4 mm locking screw - L16 mm	2
SDT2.4L18	Ø2.4 mm locking screw - L18 mm	2
SDT2.4L20	Ø2.4 mm locking screw - L20 mm	2
SDT2.4L22	Ø2.4 mm locking screw - L22 mm	1
CT2.4L14	Ø2.4 mm non-locking screw - L14 mm	1
CT2.4L16	Ø2.4 mm non-locking screw - L16 mm	1
CT2.4L18	Ø2.4 mm non-locking screw - L18 mm	1



Radial column kit – Distal radius – Size 2 – Symmetrical

KIT-WEL2S

KIT-WEL2S Radial column kit - Implants only - Distal radius - Size 2 - Symmetrical		
Ref.	Description	Qty
DETSL2	Distal radial column plate - Distal radius - Symmetrical - Size 2	1
SDT2.4L12	Ø2.4 mm locking screw - L12 mm	2
SDT2.4L14	Ø2.4 mm locking screw - L14 mm	2
SDT2.4L16	Ø2.4 mm locking screw - L16 mm	2
SDT2.4L18	Ø2.4 mm locking screw - L18 mm	2
SDT2.4L20	Ø2.4 mm locking screw - L20 mm	2
SDT2.4L22	Ø2.4 mm locking screw - L22 mm	1
CT2.4L14	Ø2.4 mm non-locking screw - L14 mm	1
CT2.4L16	Ø2.4 mm non-locking screw - L16 mm	1
CT2.4L18	Ø2.4 mm non-locking screw - L18 mm	1



Dorso-lateral kit – Distal radius – Size 2 – Right or Left

KIT-WQ2x

KIT-WQ2G OR KIT-WQ2D* Dorso-lateral kit - Implants only - Distal radius - Size 2 - Left or Right		
Ref.	Description	Qty
DTxQ2	Dorso-lateral plate - Distal radius - Size 2 - Left or Right	1
SDT2.4L12	Ø2.4 mm locking screw - L12 mm	2
SDT2.4L14	Ø2.4 mm locking screw - L14 mm	2
SDT2.4L16	Ø2.4 mm locking screw - L16 mm	2
SDT2.4L18	Ø2.4 mm locking screw - L18 mm	2
SDT2.4L20	Ø2.4 mm locking screw - L20 mm	2
SDT2.4L22	Ø2.4 mm locking screw - L22 mm	1
CT2.4L14	Ø2.4 mm non-locking screw - L14 mm	1
CT2.4L16	Ø2.4 mm non-locking screw - L16 mm	1
CT2.4L18	Ø2.4 mm non-locking screw - L18 mm	1


* Each kit includes only one plate (left, right or symmetrical)
NB : additional sterile screws can be ordered separately (S-Box)

Please note that the «x» in the reference represents the side of the implant and has to be replaced by «G» for left kits or «D» for right implants/kits.

Initial R Xpert™ 2.4 Fragment Specific - Kit references



Distal ulna kit – Size 1
Right or Left



Distal ulna kit – Size 2
Right or Left

KIT-WHE1G OR KIT-WHE1D* Distal ulna kit - Implants only - Size 1 - Left or Right		
Ref.	Description	Qty
HTxE1	Distal ulna plate - Size 1 - Left or right	1
SDT2.4L12	Ø2.4 mm locking screw - L12 mm	2
SDT2.4L14	Ø2.4 mm locking screw - L14 mm	2
SDT2.4L16	Ø2.4 mm locking screw - L16 mm	2
SDT2.4L18	Ø2.4 mm locking screw - L18 mm	2
SDT2.4L20	Ø2.4 mm locking screw - L20 mm	2
SDT2.4L22	Ø2.4 mm locking screw - L22 mm	1
CT2.4L14	Ø2.4 mm non-locking screw - L14 mm	1
CT2.4L16	Ø2.4 mm non-locking screw - L16 mm	1
CT2.4L18	Ø2.4 mm non-locking screw - L18 mm	1

KIT-WHE2G OR KIT-WHE2D* Distal ulna kit - Implants only - Size 2 - Left or Right		
Ref.	Description	Qty
HTxE2	Distal ulna plate - Size 2 - Left or right	1
SDT2.4L12	Ø2.4 mm locking screw - L12 mm	2
SDT2.4L14	Ø2.4 mm locking screw - L14 mm	2
SDT2.4L16	Ø2.4 mm locking screw - L16 mm	2
SDT2.4L18	Ø2.4 mm locking screw - L18 mm	2
SDT2.4L20	Ø2.4 mm locking screw - L20 mm	2
SDT2.4L22	Ø2.4 mm locking screw - L22 mm	1
CT2.4L14	Ø2.4 mm non-locking screw - L14 mm	1
CT2.4L16	Ø2.4 mm non-locking screw - L16 mm	1
CT2.4L18	Ø2.4 mm non-locking screw - L18 mm	1

* Each kit includes only one plate (left or right)
NB: additional sterile screws can be ordered separately (S-Box)

Please note that the «x» in the reference represents the side of the implant and has to be replaced by «G» for left kits or «D» for right implants/kits.

Removal kit Sterile instruments

REMOVAL KIT		
Ref.	Description	Content
KIT-REMOVE-2	Removal kit for T8 hexalobe	T8 prehensor screwdriver



Supplemental instrumentation kits (KIT-REMOVE-2)

INITIAL R™ XPERT - INSTRUMENTATION	
Ref.	Description
KIT-INSTRUM-1	Instrumentation kit for distal radius & distal ulna implants

INITIAL R™ XPERT - INSTRUMENTATION CONTENT	
Description	Qty
Pins - Ø1.4 L 120 mm	4
T8 prehensor screwdriver	1
Length gauge	1
Ø1.8 mm quick coupling drill bit - L 125 mm	1
Ø1.8 mm threaded guide gauge	1
Non threaded polyaxial drill guide	1



KIT-INSTRUM-1

This information is intended to demonstrate the Newclip Technics portfolio of medical devices. Always refer to the package insert, product label and/or user instructions before using any Newclip Technics product. These products must be handled and/or implanted by trained and qualified staff who have read the instructions before use. A surgeon must always rely on her or his own professional clinical judgement when deciding whether to use a particular product when treating a particular patient. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Newclip Technics representative if you have questions about the availability of Newclip Technics products in your area.



NEWCLIP TECHNICS (HQ)
45 rue des Garottières
44115 Haute Goulaine, France
+33 (0)2 28 21 23 25
orders@newcliptechnics.com
www.newcliptechnics.com

NEWCLIP TECHNICS USA
Newclip USA
340 Tesconi Circle Suite A
Santa Rosa CA 95401, USA
+1 707 230 5078
customerservice@newclipusa.com
www.newcliptechnics.com

NEWCLIP TECHNICS IBERIA
Calle Frederic Mompou 4b
Sant Just Desvern
08960 Barcelona, España
+34 938 299 526
contact@newclipiberia.com
www.newcliptechnics.com

NEWCLIP TECHNICS GERMANY
Newclip GmbH
Pröllstraße 11
D-86157 Augsburg, Deutschland
+49 (0)821 650 749 40
info@newclipgmbh.com
www.newcliptechnics.com

NEWCLIP TECHNICS JAPAN
Newclip Technics Japan K.K.
KKK Bldg. 502, 3-18-1 Asakusabashi
Taito-Ku, Tokyo, 111-0053, Japan
+81 (0)3 58 25 49 81
Fax: +81 (0)3 58 25 49 86
www.newcliptechnics.com

NEWCLIP TECHNICS AUSTRALIA
Newclip Australia
3B/11 Donkin Street
West End 4101, Australia
+61 (0)2 81 886 110
solutions@newclipaustralia.com
www.newcliptechnics.com

NEWCLIP TECHNICS BELGIUM
Derbystraat 43,
9051 Sint-Denijs-Westrem
Phone +32 (0)9 2311331
contact@newclipbelgium.com
www.newcliptechnics.com