

Initial S™ EVO Jones Fracture : Single-use sterile kit for Jones fracture screws


Initial S™ EVO Jones Fracture - Kit content

Intended purpose:

The Initial S™ single use sterile kits consist of a set of instruments or a set of instruments and implants (washers) to be used only in trauma or orthopaedic surgery, when Newclip Technics screws for osteosynthesis are fitted. Implants of the Stand-Alone Screws range are intended for fractures fixation, osteotomies and arthrodeses of bones in adults, appropriate for the size of the device.

Contraindications :

- Pregnancy.
- Acute or chronic local or systemic infections.
- Allergy to one of the materials used or sensitivity to foreign bodies.



Taps - Available separately



ANC1567-ST & ANC1568-ST



NEWCLIP-TECHNICS



Initial S™ EVO Jones Fracture - References

Initial S™ EVO Jones Fracture kit



INITIAL S™ EVO JONES FRACTURE KIT - INSTRUMENTATION CONTENT

Ref.	Description	Content	Qty
KIT-SCM4.5-5.5F	Instrumentation kit for Ø4.5 mm and Ø5.5 mm compressive screws	3 in 1 instrument for Ø4.5 mm & Ø5.5 mm screws	1
		5.8 mm single use handle	1
		Ø2.9 mm quick coupling drill bit - cannula Ø1.7 mm - L 170 mm	1
		Double guide Ø1.7 mm / Ø5.2 mm	1
		Pin Ø1.6 L150 mm	3

Implants

Ø4.5 mm & Ø5.5 mm Jones fracture screws available separately in sterile version (S-box).

Ø4.5 MM COMPRESSIVE SCREWS

Ref.	Description
MT4.5L35D-ST	Ø4.5 mm compressive screw - Partially threaded - L 35 mm - STERILE
MT4.5L40D-ST	Ø4.5 mm compressive screw - Partially threaded - L 40 mm - STERILE
MT4.5L45D-ST	Ø4.5 mm compressive screw - Partially threaded - L 45 mm - STERILE
MT4.5L50D-ST	Ø4.5 mm compressive screw - Partially threaded - L 50 mm - STERILE
MT4.5L55D-ST	Ø4.5 mm compressive screw - Partially threaded - L 55 mm - STERILE
MT4.5L60D-ST	Ø4.5 mm compressive screw - Partially threaded - L 60 mm - STERILE
MT4.5L65D-ST	Ø4.5 mm compressive screw - Partially threaded - L 65 mm - STERILE

Ø5.5 MM COMPRESSIVE SCREWS

Ref.	Description
MT5.5L35D-ST	Ø5.5 mm compressive screw - Partially threaded - L 35 mm - STERILE
MT5.5L40D-ST	Ø5.5 mm compressive screw - Partially threaded - L 40 mm - STERILE
MT5.5L45D-ST	Ø5.5 mm compressive screw - Partially threaded - L 45 mm - STERILE
MT5.5L50D-ST	Ø5.5 mm compressive screw - Partially threaded - L 50 mm - STERILE
MT5.5L55D-ST	Ø5.5 mm compressive screw - Partially threaded - L 55 mm - STERILE
MT5.5L60D-ST	Ø5.5 mm compressive screw - Partially threaded - L 60 mm - STERILE
MT5.5L65D-ST	Ø5.5 mm compressive screw - Partially threaded - L 65 mm - STERILE

This information is intended to demonstrate the Newclip Technics portfolio of medical devices. Always refer to the package insert, product label and/or user instructions before using any Newclip Technics product. These products must be handled and/or implanted by trained and qualified staff who have read the instructions before use. A surgeon must always rely on her or his own professional clinical judgement when deciding whether to use a particular product when treating a particular patient. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Newclip Technics representative if you have questions about the availability of Newclip Technics products in your area.

Taps

Sterile instruments available separately.

ASSOCIATED TAPS

Ref.	Description
ANC1567-ST	Tap for Ø4.5 mm screws - STERILE
ANC1568-ST	Tap for Ø5.5 mm screws - STERILE

CORRESPONDING TABLE

The taps allow to select the appropriate diameter of implant.



Removal kit

Sterile instruments available separately.

REMOVAL KIT

Ref.	Description	Content	Qty
KIT-REMOVE-8	Removal kit for T20 hexalobe	T20 screwdriver	1



NEWCLIP TECHNICS (HQ)

45 rue des Garotières
44115 Haute-Goulaine, France
+33 (0)2 28 21 23 25
orders@newcliptechnics.com
www.newcliptechnics.com

NEWCLIP TECHNICS IBERIA

Calle Frederic Mompou, 4b
Sant just Desvern
08960 Barcelona, España
+34 937 503 421
contact@newclipiberia.com
www.newcliptechnics.com

NEWCLIP TECHNICS USA

Newclip USA
340 Tesconi Circle, Suite A
Santa Rosa CA95401, USA
+1 707 230 5078
customerservice@newclipusa.com
www.newcliptechnics.com

NEWCLIP TECHNICS GERMANY

Newclip GmbH
Pröllstraße 11, D-86157 Augsburg,
Deutschland
+49 (0)821 650 749 40
info@newclipgmbh.com
www.newcliptechnics.com

NEWCLIP TECHNICS JAPAN

Newclip Technics Japan K.K.
KKK Bldg. 502, 3-18-1 Asakusabashi
Taito-Ku, Tokyo, 111-0053, Japan
+81 (0)3 58 25 49 81
www.newcliptechnics.com

NEWCLIP TECHNICS AUSTRALIA

Newclip Australia
3B/11 Donkin Street
West End 4101, Australia
+61 (0)2 81 886 110
solutions@newclipaustralia.com
www.newcliptechnics.com

NEWCLIP BELGIUM

Derbystraat 43,
9051 Sint-Denijs-Westrem
Phone +32 (0)9 2311331
contact@newclipbelgium.com
www.newcliptechnics.com