

INAXIS.



SET DESCRIPTION



INAXIS.

Intended purpose:

The implants of the Stand-Alone Screws range are intended for fractures fixation, osteotomies and arthrodeses of bones in adults, appropriate for the size of the device.

Contraindications:

- Pregnancy.
- Acute or chronic local or systemic infections.
- Allergy to one of the materials used or sensitivity to foreign bodies.

3 MAJOR CONFIGURATIONS FOR THE INAXIS SET

1 • COMPLETE SET



- Ø2.0 instrumentation and screws
- Ø3.0 and Ø4.0 mm instrumentation
- Ø3.0 and Ø4.0 mm screws
- Specific instrumentation

Pages 3 & 4

2 • 2 OPTIONS SET (1/2 BASE)



Choice of 2 modules between 4 options:

- Ø2.0 instrumentation and screws*
- Ø3.0 and Ø4.0 mm instrumentation
- Ø3.0 and Ø4.0 mm screws*
- Specific instrumentation

Page 5

3 • SPECIFIC INSTRUMENTATION SET (1/4 BASE)



- Specific instrumentation


Page 6

*Implants are also available in sterile


1. Complete set (1/2)



B1: Ø2.0mm instrumentation & screws Base
(ANC1572/B1 + ANC1572/R1)




B2: Ø3.0 & Ø4.0 mm instrumentation
(ANC1572/B2)




1. Complete set (2/2)



B3: Specific instrumentation (ANC1572/B3)



B4: Ø3.0 & Ø4.0 mm screws (ANC1572/B4 + ANC1572/R2 + ANC1572/R3)




2.2 options set


Choice of 2 modules between the following four:




B1: Ø2.0mm instrumentation & screws Base
(ANC1572/B1 + ANC1572/R1)




B2: Ø3.0 & Ø4.0 mm instrumentation
(ANC1572/B2)



B3: Specific instrumentation
(ANC1572/B3)




B4: Ø3.0 & Ø4.0 mm screws
(ANC1572/B4 + ANC1572/R2 + ANC1572/R3)



3. Specific instrumentation set



B3: Specific instrumentation (ANC1572/B3)



Implant & trays references not presented above.

Neutralization extremity screws

Ref.	Description
H0.9GLT2.0L12D to H0.9GLT2.0L40D	Ø2.0 mm neutralization extremity screw* - cannula Ø0.9 - full thread - L12 to 40 mm - 2mm increments
H1.3GLT3.0L25D to H1.3GLT3.0L70D	Ø3.0 mm neutralization extremity screw** - cannula Ø1.3 - full thread - L25 to 70 mm - 5mm increments
H1.3GLT4.0L25D to H1.3GLT4.0L70D	Ø4.0 mm neutralization extremity screw*** - cannula Ø1.3 - full thread - L25 to 70 mm - 5mm increments



*Not anodized
**Gold anodized
***Pink anodized

Complete set

Ref.	Description
ANC1434/B1	Set - 1 level - Size 1 - Base
ANC1573/C	InAxis set - Size 1 – Lid

2 options set

Ref.	Description
ANC1434/B4	Set - 1 level - Size 1/2 - Base
ANC1572/C2	InAxis set - Size 1/2 - Lid

1 option set

Ref.	Description
ANC1572/C1	InAxis set - Size 1/4 - Lid



This information is intended to demonstrate the Newclip Technics portfolio of medical devices. Always refer to the package insert, product label and/or user instructions including cleaning and sterilization before using any Newclip Technics product. These products must be handled and/or implanted by trained and qualified staff who have read the instructions before use. A surgeon must always rely on her or his own professional clinical judgement when deciding whether to use a particular product when treating a particular patient. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Newclip Technics representative if you have questions about the availability of Newclip Technics products in your area.

Manufacturer: Newclip Technics - Compo Kit EN - Inaxis - Ed.1 - 07/2026 - Medical devices: class IIb - CE1639 SGS BE - Read labelling and instructions before the use of Newclip Technics medical devices. These products must be handled and/or implanted by trained and qualified staff who have read the instructions before use. Non-contractual pictures.
Newclip Technics - 45 rue des Garottières - 44115 Haute Goulaine, France. Our subsidiaries: Newclip USA - Newclip Australia - Newclip GMBH - Newclip Japan - Newclip Iberia - Newclip Belgium - Newclip Italia - Newclip UK.

newcliptechnics.com

