

PERIACTIV

Indications : The implants of the PERIACTIV are intended for osteosynthesis of femoral periprosthetic fractures in adults.

Contra-indications :

- Serious vascular deterioration, bone devitalization.
- Pregnancy.
- Acute or chronic local or systemic infections.
- Lack of musculo-cutaneous cover, severe vascular deficiency affecting the concerned area.
- Insufficient bone *Altération osseuse ne permettant pas une tenue correcte des vis dans l'os.*
- Muscular deficit, neurological deficiency or behavioral disorders, which could submit th implant to abnormal mechanical strains.
- Allergy to one of the materials used or sensitivity to foreign bodies.
- Serious problems of non-compliance, mental or neurological disorders, failure to follow post-operative care recommendations.
- Unstable physical and/or mental condition.

DIAPHYSEAL PLATE

- Symmetrical plate, available in 2 sizes




DISTAL SUPRA CONDYLAR PLATE

- Asymmetrical, available in 3 sizes




TECHNICAL FEATURES

- **An anatomical plate**
The plate design follows the anatomic features of the femur, delineated notably by Pr Lefèvre (Brest).
- **Optimized fatigue resistance**
The plate is made of titanium alloy and its specific design optimizes fatigue resistance.
- **Polyaxial locking fixation system (DTS)**.....
- Possible orientation of the screw before locking
- Self-tapping screws
- Reinforced core diameter
- **A conservative treatment**
The Newclip plate is anatomically contoured and is supplied with either convergent or divergent adjustable screws.
- **Optimized fatigue resistance**
The plate is made of titanium alloy and its specific design optimizes fatigue resistance.
- **A stable plating system**
The dimensions of the plate are designed for long term resistance to the mechanical stresses the femur is subjected to.



IMPLANTS REFERENCES



DIAPHYSEAL PLATES	
Ref.	Designation
QTSM1	Periprosthetic midshaft femoral plate - size 1 - L 238 mm
QTSM3	Periprosthetic midshaft femoral plate - size 3 - L 314 mm

DISTAL PLATES	
Ref.	Designation
QTDE1	Periprosthetic distal femoral plate - Right - size 1 - L163 mm
QTGE1	Periprosthetic distal femoral plate - Left - size 1 - L163 mm
QTDE2	Periprosthetic distal femoral plate - Right - size 2 - L239 mm
QTGE2	Periprosthetic distal femoral plate - Left - size 2 - L 239 mm
QTDE3	Periprosthetic distal femoral plate - Right - size 3 - L315 mm
QTGE3	Periprosthetic distal femoral plate - Left - size 3 - L315 mm

DTS SELF-TAPPING SCREW Ø4.5 MM*	
Ref.	Designation
ST4.5L15	DTS® self-tapping screw Ø4.5 mm - L 15 mm
ST4.5L18	DTS® self-tapping screw Ø4.5 mm - L 18 mm
ST4.5L21	DTS® self-tapping screw Ø4.5 mm - L 21 mm
ST4.5L24	DTS® self-tapping screw Ø4.5 mm - L 24 mm
ST4.5L27	DTS® self-tapping screw Ø4.5 mm - L 27 mm
ST4.5L30	DTS® self-tapping screw Ø4.5 mm - L 30 mm
ST4.5L35	DTS® self-tapping screw Ø4.5 mm - L 35 mm
ST4.5L40	DTS® self-tapping screw Ø4.5 mm - L 40 mm
ST4.5L45	DTS® self-tapping screw Ø4.5 mm - L 45 mm
ST4.5L50	DTS® self-tapping screw Ø4.5 mm - L 50 mm
ST4.5L55	DTS® self-tapping screw Ø4.5 mm - L 55 mm
ST4.5L60	DTS® self-tapping screw Ø4.5 mm - L 60 mm
ST4.5L65	DTS® self-tapping screw Ø4.5 mm - L 65 mm
ST4.5L70	DTS® self-tapping screw Ø4.5 mm - L 70 mm
ST4.5L75	DTS® self-tapping screw Ø4.5 mm - L 75 mm

*Blue anodized



DIAPHYSEAL PLATES	
Ref.	Designation
QTSM1-ST	Periprosthetic midshaft femoral plate - size 1 - L 238 mm - STERILE
QTSM3-ST	Periprosthetic midshaft femoral plate - size 3 - L 314 mm - STERILE

DISTAL PLATES	
Ref.	Designation
QTDE1-ST	Periprosthetic distal femoral plate - Right - size 1 - L163 mm - STERILE
QTGE1-ST	Periprosthetic distal femoral plate - Left - size 1 - L163 mm - STERILE
QTDE2-ST	Periprosthetic distal femoral plate - Right - size 2 - L239 mm - STERILE
QTGE2-ST	Periprosthetic distal femoral plate - Left - size 2 - L 239 mm - STERILE
QTDE3-ST	Periprosthetic distal femoral plate - Right - size 3 - L315 mm - STERILE
QTGE3-ST	Periprosthetic distal femoral plate - Left - size 3 - L315 mm - STERILE

DTS SELF-TAPPING SCREW Ø4.5 MM*	
Ref.	Designation
ST4.5L15-ST	DTS® self-tapping screw Ø4.5 mm - L 15 mm - STERILE
ST4.5L18-ST	DTS® self-tapping screw Ø4.5 mm - L 18 mm - STERILE
ST4.5L21-ST	DTS® self-tapping screw Ø4.5 mm - L 21 mm - STERILE
ST4.5L24-ST	DTS® self-tapping screw Ø4.5 mm - L 24 mm - STERILE
ST4.5L27-ST	DTS® self-tapping screw Ø4.5 mm - L 27 mm - STERILE
ST4.5L30-ST	DTS® self-tapping screw Ø4.5 mm - L 30 mm - STERILE
ST4.5L35-ST	DTS® self-tapping screw Ø4.5 mm - L 35 mm - STERILE
ST4.5L40-ST	DTS® self-tapping screw Ø4.5 mm - L 40 mm - STERILE
ST4.5L45-ST	DTS® self-tapping screw Ø4.5 mm - L 45 mm - STERILE
ST4.5L50-ST	DTS® self-tapping screw Ø4.5 mm - L 50 mm - STERILE
ST4.5L55-ST	DTS® self-tapping screw Ø4.5 mm - L 55 mm - STERILE
ST4.5L60-ST	DTS® self-tapping screw Ø4.5 mm - L 60 mm - STERILE
ST4.5L65-ST	DTS® self-tapping screw Ø4.5 mm - L 65 mm - STERILE
ST4.5L70-ST	DTS® self-tapping screw Ø4.5 mm - L 70 mm - STERILE
ST4.5L75-ST	DTS® self-tapping screw Ø4.5 mm - L 75 mm - STERILE

INSTRUMENTS REFERENCES

INSTRUMENTS

Ref.	Designation
ANC119-SK	3.0 mm hexagonal screwdriver with US quick coupling system
ANC120-US	Ø4.2 mm reamer with US quick coupling system
ANC210	Length gauge for Ø4.5 mm screws
ANC211	Quick coupling drill bit Ø4.0 mm
ANC212	Ø4 mm DTS® Trauma drill guide
ANC312	Quick coupling 3.0 mm hexagonal screwdriver
ANC352	Ø6.0 mm US quick coupling handle
ANC389	QTSMx template of plate
ANC749	QTxE template of plate

The information presented in this brochure is intended to demonstrate a NEWCLIP TECHNICS product. Always refer to the package insert, product label and/or user instructions before using any NEWCLIP TECHNICS product. Surgeons must always rely on their own clinical judgment when deciding which products and techniques to use with their patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your NEWCLIP TECHNICS representative if you have questions about the availability of NEWCLIP TECHNICS products in your area.



NEWCLIP TECHNICS
PA de la Lande Saint Martin - 45 rue des Garottières
44115 Haute Goulaine (France)
Tél. : +33 (0)2 28 21 37 12 - Fax : +33 (0)2 40 63 68 37
orders@newcliptechnics.com - www.newcliptechnics.com

NEWCLIP TECHNICS AUSTRALIA
3B/11 Donkin Street
West End 4101, Australia
Phone : + 61 (0)2 81 886 110
solutions@newclipaustralia.com - www.newcliptechnics.com